



The PleuraFlow® Active Clearance Technology® System Configuration Guide


The images below represent some common configurations of PleuraFlow and conventional chest tubes with drainage tubing and canister.




One PleuraFlow and one conventional chest tube draining into one canister




Two PleuraFlow chest tubes draining into one canister




Two PleuraFlow chest tubes and one conventional chest tube draining into two canisters



Two PleuraFlow chest tubes and one conventional chest tube draining into one canister



One PleuraFlow and two conventional chest tubes draining into one canister



One PleuraFlow chest tube draining into one canister